

## Q-521 Q-Motion® Miniature Linear Stage

Smallest linear stage with position control, high resolution and attractive price



- Only 21 mm wide and 10 mm high
- Direct position measurement with integrated incremental encoder, up to 1 nm resolution (optional)
- Up to 1 nm encoder resolution
- Up to 2 nm minimum incremental motion
- Set-up of multi-axis systems with adapter plate or bracket (optionally available)
- Velocity 10 mm/s
- Suitable for vacuum to  $10^{-6}$  hPa, versions to  $10^{-9}$  hPa available

### Precision-class micropositioning stage

Q-Motion® stages are distinguished by their extremely small design and high position resolution in the nanometer range. The piezomotor drive principle and the electrical operation are cost-efficient and can be customized.

### PIShift piezo inertia drives

Self-locking when at rest, therefore no heat generation and no servo jitter. Velocity to 10 mm/s. 1 N push/pull force.

### Direct-measuring principle

Versions with noncontact optical linear encoder available. Resolution 4 nm or 1 nm, depending on the version. Versions with encoder feature a reference point switch.

### Vacuum and nonmagnetic environments

All Q-Motion® stages are suitable for operation in high vacuum to  $10^{-6}$  hPa. Furthermore, ultrahigh vacuum variants for  $10^{-9}$  hPa are also available. Nonmagnetic versions are available on request.

### Fields of application

Industry and research, metrology, microscopy, micromanipulation, biotechnology and automation.

## Compatible controllers / drivers / amplifiers

E-873 Q-Motion® servo controller  
 E-870 PIShift drive electronics  
 E-871 networkable servo controller for piezomotors  
 E-873.3QTU Q-Motion® servo controller

## Compatible accessories

Q-101.AP1 adapter plate  
 Q-121.xxx adapter bracket and adapter plate

## Related products

Q-614 Q-Motion® miniature rotation stage  
 Q-622 Q-Motion® miniature rotation stage  
 Q-632 Q-Motion® rotation stage  
 Q-545 Q-Motion® precision linear stage  
 Q-522 Q-Motion® miniature linear stage


## Specifications

	Q-521.130	Q-521.140 / Q-521.14U	Q-521.230	Q-521.240 / Q-521.24U	Q-521.330	Q-521.340 / Q-521.34U	Q-521.x00 / Q-521.x0U	Unit
<b>Motion and positioning</b>	12 mm travel range, resolution 4 nm	12 mm travel range, resolution 1 nm, UHV version Q-521.14U	22 mm travel range, resolution 4 nm	22 mm travel range, resolution 1 nm, UHV version Q-521.24U	32 mm travel range, resolution 4 nm	32 mm travel range, resolution 1 nm, UHV version Q-521.34U	12 mm to 32 mm travel range, open-loop, UHV versions Q-521.x0U	
Active axis	X	X	X	X	X	X	X	
Travel range	12	12	22	22	32	32	12 to 32	mm
Integrated sensor	Linear encoder	Linear encoder	Linear encoder	Linear encoder	Linear encoder	Linear encoder	-	
Sensor resolution	4	1	4	1	4	1	-	nm
Min. incremental motion	8	2	8	2	8	2	-	nm
Unidirectional repeatability over entire travel range	25	25	25	25	30	30	-	nm
Bidirectional repeatability over entire travel range	40	40	40	40	50	50	-	nm
Unidirectional repeatability over 100 µm travel range	12	12	12	12	12	12	-	
Bidirectional repeatability over 100 µm travel range	24	24	24	24	24	24	-	nm
Pitch / yaw over entire travel range	100	100	100	100	100	100	100	µrad
Pitch / yaw over 100 µm travel range	1	1	1	1	1	1	1	µrad
Maximum velocity*	10	10	10	10	10	10	10	mm/s


	Q-521.130	Q-521.140 / Q-521.14U	Q-521.230	Q-521.240 / Q-521.24U	Q-521.330	Q-521.340 / Q-521.34U	Q-521.x00 / Q-521.x0U	Unit
<b>Mechanical properties</b>								
Load capacity	10	10	10	10	10	10	10	N
Length	30	30	32.2	32.2	42.2	42.2	30 to 42.2	mm
Width	21	21	21	21	21	21	21	mm
Height	10	10	10	10	10	10	10	mm
<b>Drive properties</b>								
Motor type	Piezoelectric inertia drive	Piezoelectric inertia drive	Piezoelectric inertia drive	Piezoelectric inertia drive	Piezoelectric inertia drive	Piezoelectric inertia drive	Piezoelectric inertia drive	
Push / pull force	1	1	1	1	1	1	1	N
Holding force, de-energized	1.3	1.3	1.3	1.3	1.3	1.3	1.3	N
<b>Miscellaneous</b>								
Operating temperature range	0 to 40	0 to 40	0 to 40	0 to 40	0 to 40	0 to 40	0 to 40	°C
Material	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	
Cable length	1	1	1	1	1	1	1	m
Connector	Sub-D	Sub-D	Sub-D	Sub-D	Sub-D	Sub-D	Sub-D	
Recommended controller	E-871, E-873	E-871, E-873	E-871, E-873	E-871, E-873	E-871, E-873	E-871, E-873	E-870	

\* Typical velocity at a control frequency of 20 kHz


## Drawings and images



Q-521.1xx, Dimensions in mm. Note that the decimal places are separated by a comma in the drawings.



Q-521.2xx, Dimensions in mm. Note that the decimal places are separated by a comma in the drawings.



Q-521.3xx, Dimensions in mm. Note that the decimal places are separated by a comma in the drawings.

## Ordering information

### Q-521.100

Q-Motion® miniature linear positioning stage, 12 mm travel range, without position sensor for open-loop operation, 1 N push/pull force, dimensions 21 × 30 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.130

Q-Motion® miniature linear positioning stage, 12 mm travel range, linear encoder, 4 nm resolution, 1 N push/pull force, dimensions 21 × 30 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.140

Q-Motion® miniature linear positioning stage, 12 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 30 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.14U

Q-Motion® miniature linear positioning stage, 12 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 30 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-9</sup> hPa

### Q-521.200

Q-Motion® miniature linear positioning stage, 22 mm travel range, without position sensor for open-loop operation, 1 N push/pull force, dimensions 21 × 32 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.230

Q-Motion® miniature linear positioning stage, 22 mm travel range, linear encoder, 4 nm resolution, 1 N push/pull force, dimensions 21 × 32 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.240

Q-Motion® miniature linear positioning stage, 22 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 32 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.24U

Q-Motion® miniature linear positioning stage, 22 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 32 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-9</sup> hPa

### Q-521.300

Q-Motion® miniature linear positioning stage, 32 mm travel range, without position sensor for open-loop operation, 1 N push/pull force, dimensions 21 × 42 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

### Q-521.330

Q-Motion® miniature linear positioning stage, 32 mm travel range, linear encoder, 4 nm resolution, 1 N push/pull force, dimensions 21 × 42 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

**Q-521.340**

Q-Motion® miniature linear positioning stage, 32 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 42 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-6</sup> hPa

**Q-521.34U**

Q-Motion® miniature linear positioning stage, 32 mm travel range, linear encoder, 1 nm resolution, 1 N push/pull force, dimensions 21 × 42 × 10 mm (W × L × H), piezoelectric inertia drive, vacuum-compatible to 10<sup>-9</sup> hPa

Ask about custom designs!