

**1. 01 T30X-6C**      6 mm Travel, Miniature Linear Translation Stage, Center Driving

**Specifications**

Travel	6 mm
Platform Dimensions	30 mm x 30 mm
Actuator	Micrometer Drive in Center
Travel Guide	Cross Roller Bearing
Material	Aluminum
Finish	Black Anodized
Max. Load	1.5 kg
Weight	0.06 kg



Minimum Readout	0.01 mm
Accuracy	0.01 mm
Degree of Parallelization	0.02 mm
Part Number	T30X-6C

**Features**

- Heavy-Duty Cross Roller Bearing
- Micrometer Drive
- M2 Mounting Holes

**Dimension**

