


HPG (HIGH POWER GREEN) SERIES LASER



Model Name	HPG-2000	HPG-3000	HPG-4000	HPG-5000	HPG-6000
Wavelength	532nm				
Output Power	2W	3W	4W	5W	6W
Spatial Mode	TEM ₀₀ , M ² < 1.2				
Power Stability	<+/- 2 % over 1 hour				
Noise	3% Pk/Pk				
Beam Diameter	<2mm at output aperture				
Beam Divergence	< 1.5mRad				
Polarization	>100:1, Horizontal				
Laser Head Size	132mm x 110mm x 376mm				
Power supply	2U high 19 inch rack				
Control	Local & RS232				
Cooling	Air cooled				
Power requirement	100 - 240 VAC, 50 / 60 Hz Max 1000W				



Elforlight Ltd, 4b Brunel Close, Daventry, Northants NN11 8RB, UK

Tel +44(0)1327 300069 | Fax +44(0)1327 300076

Web: www.elforlight.com | Email: info@elforlight.com