LASER SYSTEMS

Gigashot™FT

FEATURES AND BENEFITS

PULSED Nd: YAG DPSS LASER SYSTEM


GIGASH®

- Smooth flat top profile
 - 320 mJ @ 1064 nm, < 10 nsec, 100 Hz
- 532 nm / 355 nm available
 - Injection seeding option
 - Low maintenance
 - High efficiency
 - Long life diode bars
- Beam characteristics maintained over a wide adjustable operating energy range
- eDrive[™] control electronics with digital remote control


The Gigashot™ FT is a new diode-pumped solid–state (DPSS) Nd:YAG laser that delivers 320 mJ per pulse at 1064 nm at a repetition rate of 100 Hz. The laser has a Master Oscillator-Power Amplifier (MOPA) architecture with an output beam that has a 'flat top' beam profile in the near field – making it ideally suited for pumping ultrafast Ti:Sapphire amplifiers and OPCPAs.

Diode pumping delivers excellent pulse to pulse and long term stability as well as a long operating lifetime. The laser is supplied with a 2 year/10,000 operating hour warranty. Injection seeding is optional.

Customized versions of the Gigashot $^{\text{TM}}$ laser system are also available. Please contact CEO® for more information.



Near field beam profile, 532 nm at 100 Hz




Near field beam profile, 355 nm at 100 Hz


LASER SYSTEMS

Gigashot™FT

GIGASHOT SPECIFICATIONS

arameter	Configurations			Units
Model	GS-320-QFI	GS-160-QFG	GS-120-QFU	_
aser Type	DPSS Nd:YAG	DPSS Nd:YAG	DPSS Nd:YAG	_
Vavelength	1064	532	355	nm
Repetition Rate	100	100	100	Hz
Output	> 320	> 160	> 120	mJ
Spatial Mode*	< 15	< 15	< 15	% rms
Beam Diameter @ Output Window	< 6.5	< 6.5	< 6.5	mm
Beam Divergence (Full Angle)	< 0.6	< 0.6	< 0.6	mrad
Pulse Width (FWHM)	< 10	< 10	< 10	nsec
Pulse-to-Pulse Energy Stability	< 1.0	< 1.5	< 2	% rms
Jitter**	< 1.0	< 1.0	< 1.0	ns
Output Stability Over 8 hr	< 2	< 2	< 2	% rms
Polarization	Vertical	Horizontal	Vertical	_
Electrical @ 50/60 Hz (Auto Ranging)	85 - 264	85 - 264	85 - 264	VAC
Operating Temperature (non-condensing)	18 - 30°C	18 - 30°C	18 - 30°C	°C
Dimensions	56.00 x 34.00 x 10.77	56.00 x 34.00 x 10.77	56.00 x 34.00 x 10.77	in





Dimensions in Inches [cm]

Copyright © 2015 Northrop Grumman Cutting Edge Optronics All Rights Reserved. Northrop Grumman Cutting Edge Optronics reserves the right to change product design and specifications at any time without notice. No license is granted by implication or otherwise under any patents or patent rights of Northrop Grumman Cutting Edge Optronics or others. No responsibility is assumed for the use of these products, nor for any infringement on the rights of others resulting from the use of these products. Information contained herein is believed to be reliable and accurate. Laser diode product components are intended for use in a user-devised end system. However, these products are capable of emitting Class IV radiation. Extreme care must be exercised during their operation. Only persons familiar with the appropriate safety precautions should operate a laser product. Directly viewing the laser beam or exposure to specular reflections must be avoided. Serious injury may result if any part of the body is exposed to the beam. The eye is extremely sensitive to the infrared radiation and therefore, proper eyewear must be worn at all times. Use of optical instruments with these products may increase eye hazard. Always wear proper eye protection when operating.

This Product is covered by one or more of the following Patents: 5,898,211 5,985,684 5,913,108 6,310,900 Other US and Foreign Patents Pending.



5.08 [129] BEAM HEIGHT